themelves with implantology techniques. Improper technique can cause bone loss and implant failure. MOR® Implants are intended to be used only with MOR® Implant System. Some designs use bone drills and prosthetics attachments. Implants placed at severe angles relative to existing dentition or multiple implants placed at convergent/divergent manner can result in complex restorations that may overload implants, potentially leading to implant failure. A thorough diagnostic work-up and use of surgical templates are recommended to help ensure proper angulation.

Other relative contraindications include steroid and anticoagulant treatment which may affect the surgical site, surrounding tissue, or patient's healing function. Exposure to long-term use of bisphosphonate drugs especially with chemotherapy may impact implant survival. Patient selection including consultation with the attending physician is strongly recommended prior to treatment. Excessive mobility, bone loss, or infection may indicate the implant is failing. Any implant which appears to be failing should be treated or removed as soon as possible. If removal is necessary, curette any soft tissue from the implant site and allow site to heal as though it were an traumatic extraction. Due to the metal conductivity, electrosurgery around the implants and intraoral abutment preparations without irrigation will result in tissue damage and implant failure. Patients should consult with their physician and imaging technician prior to undergoing an MRI procedure.

Adequate preoperative imaging and diagnostic evaluation are necessary to determine available bone anatomy in prospective implant sites. A thorough clinical examination and assessment of complete medical and dental history to determine whether or not a patient is a good implant candidate.

2. Communication and planning with each clinician and laboratory technician who is involved in the treatment plan of the patient.

3. Mounted diagnostic casts.
***NOTE: MOR® Implant surgery may be performed as a guided technique as well which always recommended.**
When using a digitally planned and fully guided technique, follow the specific sequence of steps indicated by the planning software in coordination with the Sterngold MOR® surgical instruments and implants.

Clinicians choosing to perform a non-guided procedure should assess the following:
a. CT Scans/Radiographs for assessing bone quality, anatomy of the ridge, maxillary sinus, nasal cavities, inferior alveolar nerve, mental foramen, natural tooth positions and other anatomic features that may affect implant placement.
b. An evaluation of mandibular bone morphology and skeletal jaw classification.
c. Determination of the type, size, quantity, and location of the implants for optimum esthetics and functional results.
d. Fabrication of a surgical template which closely duplicates the position of the teeth in the final restoration or a template into which drills are positioned for the osteotomy and final implant placement is fabricated.

IMPLANT SURGERY
1. Administer anesthesia in the appropriate manner.
2. Whether the MOR® Implants are to be used on a temporary or long-term basis, the MOR® Attachment Head of the implant will protrude through the tissue. Therefore, when placing the MOR® Implants, drilling through the tissue, a tissue punch, or a small flap may be used. To use a tissue punch, the ridge must be of adequate width to accommodate the diameter of the head of the MOR® Implants, plus a reasonable amount of additional bone. Also, adequate attached gingiva must be present. If the MOR® Implants are to be placed in the anterior or if you are also placing additional implants or grafting, it is necessary to identify the neurovascular bundles of the mental foramen in the mandible and the borders of the maxillary sinus in the maxilla. This is done by making an incision through the mucosa and attached gingiva along the crest of the ridge and deflecting a full-thickness mucoperiosteal flap both lingually (palatally) and facially to expand the operative site. Examine the geometry, quality, and quantity of available bone at the site.

When placing the 2.1 mm or the 2.4 mm MOR® Implants, the MOR® Drill is used to gently penetrate the tissue at each surgical site and to penetrate the cortical plate. In most cases no incision is necessary. The depth of the osteotomy, on average, is one-third the length of the implant. When the length of the implant is 3.0 mm, the depth of the osteotomy is three-quarters to the full length of the threaded part of the implant, depending on the density of the bone. As an option, in very dense bone, the 2.7 mm drill may be used to engage the osteotomy to a depth of approximately one-half the length of the threaded part of the implant. Drilling is performed at low speeds, approximately 1200 RPM. The site is irrigated to limit the heat transfer to the bone. The implants should be placed within 10 degrees of each other and within 10 degrees of the appropriate path of insertion of the prosthesis.

4. The implant should always be advancing into the bone while it is turning. It is very dense bone and it is possible to drill the end of the MOR® Implant to not be able to pierce the undrilled bone, but the implant may still be turning. If that occurs, STOP IMMEDIATELY. Back the implant out of the osteotomy and extend the length of

the hole.
5. When placing the MOR® Implants in less dense bone, such as that often found in the maxilla, you may wish to reduce the drill speed to less than 1000 RPM to improve tactile sense and to help to avoid oversizing the osteotomy. The use of the 2.4 mm diameter MOR® Implant, instead of the 2.1 mm, is recommended, while still using the MOR® Drill to prepare the osteotomy. A shallow pilot hole is created, basically penetrating the cortical plate. Using this technique, the osteotomy will be made in the area immediately around the implant as it is being inserted.
6. To remove the implant from the packaging and drive it into the bone:
a. Open the outer box, peel back the outer cover on the blister pack, and drop the inner blister package onto the sterile field. The outer box and outer blister pack should be opened by a non-sterile assistant.
c. Use the appropriate MOR® tool (MOR® Insertion Tool or MOR® Latch Driver – Captive) to deliver the implant from the package to the implant site and to begin hand-driving the implant into the bone. Push the driver into the head of the implant. The O-ring will hold the implant in place. When torque resistance is encountered with the tool, you may switch to the MOR® Wing Driver.

d. If the insertion cannot be completed using the MOR® Wing Driver, the handpiece driven MOR® Latch Driver – Captive is placed into a surgical handpiece and engaged onto the square of the implant head. Drive the implant at slow speed (15 RPM) with irrigation, and a torque setting no higher than 45 Ncm

e. The final turn or two of the implant placement may be done by hand, using a torque wrench. If your torque wrench accepts latch-type tools, use the MOR® Latch Driver – Captive. If your torque wrench takes 4 mm square drive tools, use the MOR® Socket Wrench to complete the seating of the implant so that the square drive area of the implant head is just supragingival. Do not exceed 45 Ncm.

PRE-OPERATORY TREATMENT PLANNING:
Pre-Operative planning should consist of all of the following:
1. A thorough clinical examination and assessment of complete medical and dental history to determine whether or not a patient is a good implant candidate.
2. Communication and planning with each clinician and laboratory technician who is involved in the treatment plan of the patient.
3. Mounted diagnostic casts.
***NOTE: MOR® Implant surgery may be performed as a guided technique as well which always recommended.**
When using a digitally planned and fully guided technique, follow the specific sequence of steps indicated by the planning software in coordination with the Sterngold MOR® surgical instruments and implants.

Clinicians choosing to perform a non-guided procedure should assess the following:
a. CT Scans/Radiographs for assessing bone quality, anatomy of the ridge, maxillary sinus, nasal cavities, inferior alveolar nerve, mental foramen, natural tooth positions and other anatomic features that may affect implant placement.
b. An evaluation of mandibular bone morphology and skeletal jaw classification.
c. Determination of the type, size, quantity, and location of the implants for optimum esthetics and functional results.
d. Fabrication of a surgical template which closely duplicates the position of the teeth in the final restoration or a template into which drills are positioned for the osteotomy and final implant placement is fabricated.

CONTRAINDICATIONS:
The following conditions would obviate the use of dental implants: inadequate bone substance or quality, intractable poor oral hygiene, acute or chronic infection, abuse of drugs or alcohol, cirrhosis, allergies to titanium, smoking, cardiopathy, neoplasia in action, systemic conditions that would impair healing, excessive occlusal parafunction, history of radiation or patient otherwise not suitable for long or complicated surgery, or the inability to construct a functional prosthesis.
▲PATIENT PRECAUTIONS:
Following the surgical procedure and implant placement, the patient must carefully follow specific post-surgical directions including avoidance of prolonged strenuous physical activity following the procedure, consuming only soft foods, pain management and follow up as prescribed by the clinician.

▲WARNINGS:
MOR® Implants should not be placed if there is an insufficient volume of alveolar bone to minimally support the implant (minimum 1mm circumferential and 2mm apical). Implants placed in the maxilla should not perforate the sinus floor membrane. Poor bone quality, poor patient oral hygiene, heavy tobacco use, uncontrolled systemic diseases (diabetes, etc.), reduced immunity, alcoholism, drug addiction, and psychological instability may contribute to lack of integration and/or subsequent implant failure. Severe bruxism, clenching, and overloading may cause bone loss, screw loosening, component fracture, and/or implant failure. Exposure to radiation and chemotherapy may impact the health of the implant. Dental implant patients should be instructed to consult with their physician prior to undergoing such treatment options.

▲PRECAUTIONS - IMPLANTS:
Surgical and restorative techniques required to place dental implants are highly specialized and complex procedures. Practitioners should attend courses of study to familiarize

the hole.

5. When placing the MOR® Implants in less dense bone, such as that often found in the maxilla, you may wish to reduce the drill speed to less than 1000 RPM to improve tactile sense and to help to avoid oversizing the osteotomy. The use of the 2.4 mm diameter MOR® Implant, instead of the 2.1 mm, is recommended, while still using the MOR® Drill to prepare the osteotomy. A shallow pilot hole is created, basically penetrating the cortical plate. Using this technique, the osteotomy will be made in the area immediately around the implant as it is being inserted.

6. To remove the implant from the packaging and drive it into the bone:
a. Open the outer box, peel back the outer cover on the blister pack, and drop the inner blister package onto the sterile field. The outer box and outer blister pack should be opened by a non-sterile assistant.
c. Use the appropriate MOR® tool (MOR® Insertion Tool or MOR® Latch Driver – Captive) to deliver the implant from the package to the implant site and to begin hand-driving the implant into the bone. Push the driver into the head of the implant. The O-ring will hold the implant in place. When torque resistance is encountered with the tool, you may switch to the MOR® Wing Driver.

d. If the insertion cannot be completed using the MOR® Wing Driver, the handpiece driven MOR® Latch Driver – Captive is placed into a surgical handpiece and engaged onto the square of the implant head. Drive the implant at slow speed (15 RPM) with irrigation, and a torque setting no higher than 45 Ncm

e. The final turn or two of the implant placement may be done by hand, using a torque wrench. If your torque wrench accepts latch-type tools, use the MOR® Latch Driver – Captive. If your torque wrench takes 4 mm square drive tools, use the MOR® Socket Wrench to complete the seating of the implant so that the square drive area of the implant head is just supragingival. Do not exceed 45 Ncm.

IMPLANT PROSTHETIC RESTORATION:
Before using the Sterngold Prosthetic Components, the clinician in charge should carefully study the indications, contraindications, recommendations, warnings and instructions, as well as all other product-specific information (technical product description, catalogue sheet, etc.) and fully comply with them. Detailed instructions over and above those contained in this instruction for use concerning the possible combinations, product-specific risks, preparatory steps, indications and contraindications, etc. can be found in the description of the prosthetic instructions for use and technical user's guide. It is also recommended to attend the appropriate user-training courses. The aforementioned documents and details of the training courses may be obtained from the appropriate representatives in the various countries. The manufacturer, the importer and the suppliers of MOR® Implant System products are not liable for complications, other negative effects or damages that might occur for reasons such as incorrect indications or surgical technique, unsuitable choice of material or handling thereof, unsuitable use or handling of the instruments, asepsis and so on.

The clinician is responsible for any such complications or other consequences. It is also the clinician's responsibility to properly instruct and inform the patient on the functions, handling and necessary care of the product and on all known product risks.

INDICATIONS/INTENDED USE – PROSTHETICS COMPONENTS:
The Sterngold Dental prosthetic components are devices used to connect Sterngold dental prostheses. They may be used to be used in single and multiple implant fixed and removable restorations.
▲WARNINGS-PROSTHETICS:
The Sterngold Prosthetic Components should not be used unless the dental implants are stable and there are no signs of infection or severe bone loss around the implant.
Based on the treatment plan, prepare the indicated restorative prosthetic component.

MOR® restorative components include the following: MOR® RC-1 and RC-2 retention caps, MOR® impression coping and MOR® MOR® A TempORIZATION cap, MOR® Implant Analog and MOR® Scan Body, MOR® Comfort Cap, and Black and Extra Strong Purple O-rings.

When making a final physical impression for the prosthetic restoration, use the MOR® Impression Coping and corresponding MOR® Implant Analogs for fabrication of master casts. When making a digital scan, the MOR® Scan Body is affixed to the implant before scanning.

Implant Loading and Overdenture Pickup
After each implant has been placed and minimum initial stabilization torque (35 Ncm) has been established and overall implant stability has been determined by clinician, MOR® RC-1 or RC-2 retention caps and/or assembly is placed on each implant ready for pick up. Ensure denture housing bases created to allow the pickup of denture over the housings when fully seated in the mouth. Place a small amount of pick-up material (i.e. Sterngold EZ-PickUp®, See

IFU for details on denture preparation) on top of RC-1/RC-2 housing and fill recess ¾ full with pick-up material. Seat the denture and allow pickup material to polymerize. Clean up and prepare denture for final delivery.

▲Prosthetic Precautions:
During prosthetic placement procedures, care should be taken to minimize tissue damage and infection. Radiographs and other diagnostic reviews are needed during prosthetic treatment planning and reconstruction to ensure proper seating of implant/overdenture housing interface. Prosthodontic procedures must consider proper stress distribution, passive adjustment and fitting of prosthesis to fixture abutments, adapting of occlusion to opposite jaw, and avoiding excessive transverse load.

▲CLEANING AND STERILIZATION:
Disinfection and sterilization procedures should conform to OSHA or local guidelines for blood borne pathogens.

STERILE SINGLE-USE COMPONENTS:
All MOR® Implants have been gamma radiation sterilized and are double blister packaged for single use only. Do not re-use or re-sterilize implants. Do not use sterile devices if the packaging providing the sterile barrier is damaged or compromised. If your manner. Reuse of a single use device that has come in contact with blood, hard or soft tissue or other body fluids may lead to patient or user injury. Possible risks associated with reuse of a single use device include, but are not limited to mechanical failure and transmission of infectious agents.

NON-STERILE RE-USABLE AND SINGLE USE COMPONENTS:
All non-sterile, re-usable tools and components must be cleaned and sterilized prior to use. These include: Implant Driver Wrench, Implant Driver Latch, Drill Extender, Insertion Tool, Wing Driver, Graduated Torque Wrench, and the Kit Box, as well as all single use, non-sterile components used in the preparation of the implant. Analog, Scan Bodies, O-Rings, TempORIZATION Caps, Retaining Caps, and Comfort Cap.

CLEANING/DISINFECTION PRIOR TO USE:
Devices are to be (disassembled, if necessary, i.e., wrenches, Kit Box, etc.) cleaned of biological matter and fluids, disinfected, and sterilized by user prior to use. Use the following guidelines for cleaning/disinfecting for non-sterile re-usable components. Rinse with cool-to-lukewarm water for two-and-one-half minutes. Immerse all instruments in an ultrasonic cleaner with a neutral enzymatic detergent disinfection solution per the manufacturer's guidelines. Sonicate for 10 minutes. Rinse with tap water for three minutes. Note: Using detergents other than directed, including alcohol or acetates, may damage the Kit Box components.

TRANSPORT:
Carry implants in sterile package to sterile field without compromising sterility.

INSPECTION AND FUNCTIONALITY:
Prior to each use, special attention must be given to each instrument to confirm that it has been cleaned, visually inspected, and functionally tested to ensure it meets performance requirements as determined by the clinician. Once it has been determined that the instrument no longer meets these functional requirements, the instrument must be replaced.

All non-sterile, single-use components must be sterilized prior to use, with the exception of non-metal components, which should only be cleaned and disinfected.

STERILIZATION:
Individual parts should be reassembled as required, packaged and placed in appropriate autoclave for moist heat sterilization or in a dry heat pouch for dry heat sterilization. The following sterilization parameters (method, time, and temperature) are required to achieve a 10-6 sterility assurance level (SAL). Local or national specifications should be followed where steam sterilization requirements are stricter or more conservative than those listed in the table. Verify the calibration of your unit to ensure recommended temperatures are not exceeded. To ensure the autoclave is performing effectively, periodic use of biologic indicators should be considered.

Parts Individually Produced	Cycle Type	Temperature	Exposure Time	Dry Time (only for kits)
Parts in a kit	Gravity (steam)	121°C 250°F	40 minutes	n/a
	Gravity (steam)	132°C 270°F	15 minutes	30 minutes
	Pre-vacuum (steam)	134°C 273°F	5 minutes	30 minutes

Warning: The use of non-sterile components may lead to infection of hard and/or soft tissue or may cause infectious disease or implant failure.

▲BREAKAGE:
Implant fractures can occur when applied loads exceed the normal functional design tolerances of the implant components. Potential causes of implant failure include, but are not limited to: implant numbers, lengths and/or diameters to adequately support a restoration, excessive cantilever length, incomplete abutment seating, abutment angles greater than 0° degrees, occlusal interferences causing excessive lateral forces, patient parafunction (e.g. bruxing, clenching), improper denture manufacture procedures, inadequate denture fit, and physical trauma.

CHANGES IN PERFORMANCE:

It is the responsibility of the clinician to instruct the patient on all appropriate contraindications, side effects, and precautions as well as the need to seek the services of a trained dental professional if there are any changes in the performance of the implant (e.g., looseness of the prosthesis, infection or exudate around the implant, pain, or any other unusual symptoms that the patient has not been told to expect).

HYGIENE & MAINTENANCE:
Long-term potential health is directly related to the maintenance of oral hygiene. Potential implant candidates should establish an adequate oral hygiene regimen prior to implant therapy. Following implant placement, the clinician should instruct the patient on proper tools and techniques to ensure long-term maintenance of the implant(s). The patient should also be instructed to maintain routinely scheduled prophylaxis and evaluation appointments.

SHELF LIFE:
The product expiration date is indicated by the hourglass symbol on the product label, followed by the year and month of expiration for all implants.

▲CAUTION:
Federal (USA) law restricts this device to sale by or on the order of a dentist (or other licensed practitioner).

HOW SUPPLIED:
Sterngold Dental implants are supplied sterile and should be handled with sterile titanium instruments. Implants are intended for single use only and should not be re-sterilized if contaminated.

STORAGE:
The expiration date is provided on the product labeling. Do not use a product that has exceeded the expiration date or has a damaged sterile barrier.

DISPOSAL:
Disposal of these components shall follow government regulations and environmental requirements.

EMC AND ELECTRICAL SAFETY:
The MOR® Implants do not require EMC and Electrical Safety Evaluation.

MR SAFETY:
The MOR® Implants have not been evaluated for safety and functionality in the MR environment. The MOR® Implants have not been tested for heating, migration, or image artifact in the MR environment. The safety of MOR® Implants in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

WARRANTY:
Sterngold guarantees its products to be free from defects in material and/or workmanship. No other warranty is expressed or implied. This warranty applies to the original purchaser only. In the unlikely event of a defect, please follow the returned goods policy outlined in the product catalog.

Symbols:
The following table describes the symbols that may be printed on the packaging label. Please refer to the packaging label for the applicable symbols related to the product.

Symbol	Used For	Symbol	Used For
	Use by Prescription Only		Part No.
	Caution, consult accompanying documents		Batch Code (lot #)
	Do Not re-use		Do not use if package is damaged
	Use by YYYY-MM-DD or YYYY-MM		Do not resterilize
	Manufacturer		Sterilized using irradiation

Manufactured and Distributed in the USA by:

Sterngold Dental, LLC
23 Frank Mossberg Drive
Attleboro, MA 02703 USA
+1.800.243.9942 / +1.508.226.5660

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