or quality, infection, poor patient oral hygiene or other co-morbid conditions. Temporary or permanent numbness (anesthesia), paresthesia or dysesthesia are possible complications of implant surgery. Loosening or fracture of implant and restoration components may occur over time. Possible trauma or injury to implant site if incorrect selection of drill is used or used incorrectly, see below under Cautions.

▲ PATIENT PRECAUTIONS:

Following the surgical procedure and implant placement, the patient must carefully follow specific post-surgical directions including avoidance of prolonged strenuous physical activity following the procedure, consuming only soft foods, pain management and follow up as prescribed by the clinician.

▲ PRECAUTIONS:

Implant surgery is a highly specialized and complex procedure and special training is required. Clinicians should attend advanced courses designed to be educated in the required techniques. Note: 100% implant success is not guaranteed. Proper treatment planning and patient selection are critical for ensuring success. Improper technique and prosthetic planning can result in occlusal overload, implant failure and loss of surrounding bone. Radiographs or other diagnostic review should be performed to determine the quality and density of bone, overall condition of the maxillary sinus, nasal cavities, inferior alveolar nerve, mental foramen, natural tooth positions and other anatomical features that may affect implant placement and negatively affect the overall case success. Dutiful and accurate communication between the surgeon, restorative dentist, and dental laboratory is also essential for success. During the implant surgical procedure, special attention to thermal and surgical trauma must be given to minimize tissue damage, infection, and excessive post-operative bleeding. Thermal trauma severely impedes implant integration to the bone. Reduced drill speeds of 1000-1500 r.p.m., the use of sharp drills, sufficient irrigation, and use of pilot drills in successively increasing sizes are essential. Careful consideration should be given to the density of the bone to avoid undue pressure to the bone upon insertion which may result in damage to the implant site which could result in pressure necrosis and implant failure. When placing the implant, never exceed 55 Ncm torque upon insertion.

Pre-Operatory Treatment Planning and Surgical Pre-Operative planning should consist of all the following:

- A thorough clinical examination and assessment of complete medical and dental history to determine whether a patient is a good implant candidate.
- Communication and planning with each clinician and laboratory technician who is involved in the treatment plan of the patient.
- Mounted diagnostic casts.

- CT Scans/Radiographs for assessing bone quality, anatomy of the residual ridge, maxillary sinus, nasal cavities, inferior alveolar nerve, mental foramen, natural tooth positions and other anatomic features that may affect implant placement.
- Communication and planning with each clinician and laboratory technician who is involved in the treatment plan of the patient.

- An evaluation of bone morphology and skeletal jaw classification.
- Determination of the type, size, quantity, and location of the implants for optimum esthetics and functional results.

- Fabrication of a surgical template based on treatment planning may be recommended which closely duplicates the position of the teeth in the final restoration or a template into which drills are positioned for the osteotomy to facilitate the final implant placement.

Refer to the TRU®/PUR® instructions for use for additional instructions on Sterngold implants and implant placement procedure.

▲ CLEANING AND STERILIZATION: Disinfection and sterilization procedures should conform to OSHA or local guidelines for blood borne pathogens.

NON-STERILE REUSABLE AND SINGLE USE COMPONENTS: All non-sterile, re-usable instruments must be cleaned and sterilized prior to use. These include Implant Drivers, Surgical Drills, Bone Taps, Pilot Drills, Lance Pilot Drills, Implant Driver Wrenches, Implant Driver Latch, Drill Extender, Ø50 Hex Driver, Ø50° Hex Tools, Square to Latch Adapter, Thumb Knobs, Graduated Torque Wrench, Torque Wrench, Torque Wrench Extension, Ratchet Wrench, Ratchet Extension, Lindemann Drill, Drill Extender, Dense Bone Drills, and Crestal Bone Drills.

INDICATIONS FOR USE: The Sterngold TRU®/PUR® Surgery Kit are to be used solely with Sterngold TRU®/PUR® implants. Patient selection should be such that there are no compromising health conditions, disease, or medications that would interfere with the complete integration of the implants into the bone. Adequate bone quality and quantity should be assessed as well as available vertical opening to allow for kit instruments to be used.

CONTRAINDICATIONS: The following conditions would obviate the use of dental implants: inadequate bone quality and volume, intractable poor oral hygiene, acute or chronic infection, allergy to medical grade titanium, abuse of drugs or alcohol, hypertension, cardiopathy, neoplasia in action, systemic conditions that would impair healing, excessive occlusal parafunction, history of radiation or patient otherwise not suitable for long or complicated surgery, or the inability to construct a functional prosthesis.

▲ WARNINGS: Ensure that each drill is used as directed by the surgical plan. Failure to do so may result in permanent damage and injury to the surgical site, nerves, or other vital anatomical structures. It is possible that an implant may fail to integrate with the surrounding bone, which would lead to loss of the implant and restoration supported by the implant. Potential causes for failure include inadequate surgical technique¹ not following the directions for use, poor case planning and/or patient selection which may include lack of bone quantity

assurance level (SAL). Local or national specifications should be followed where steam sterilization requirements are stricter or more conservative than those listed in the table. Verify the calibration of your unit to ensure recommended temperatures are not exceeded. To ensure the autoclave is performing effectively, periodic use of biologic indicators should be considered.

| Pieces Individually Produced | Cycle Type | Temperature | Exposure Time | Dry Time (only for kits) |
|------------------------------|--------------------|-------------|---------------|--------------------------|
| | Gravity (steam) | 121°C 250°F | 40 minutes | n/a |
| Pairs in a kit | Gravity (steam) | 132°C 270°F | 15 minutes | 30 minutes |
| | Pre-vacumm (steam) | 134°C 273°F | 5 minutes | 30 minutes |

Warning: The use of non-sterile components may lead to infection of hard and/or soft tissue or may cause infectious disease or implant failure.

▲ BREAKAGE: TRU®/PUR® Surgical instruments should only be used for their proper intended use to avoid damage or breakage.

CHANGES IN PERFORMANCE:

It is the responsibility of the clinician to instruct the patient on all appropriate contraindications, side effects, and precautions as well as the need to seek the services of a trained dental professional if there are any changes in the performance of the implant (e.g., looseness of the prosthesis, infection or exudate around the implant, pain, or any other unusual symptoms that the patient has not been told to expect).

SHELF LIFE: The product expiration date is indicated by the hourglass symbol on the product label, followed by the year and month of expiration.

▲ CAUTION: Federal (USA) law restricts this device to the sale by or on the order of a licensed clinician).

HOW SUPPLIED: Sterngold TRU®/PUR® Surgical Kit and instruments are supplied non-sterile and should sterilized prior to use according to the above directive.

SOFTWARE: These devices do not contain or utilize software. .

STORAGE: Place devices in a dry place to prevent damage and/or deterioration.



DISPOSAL: Disposal of these components shall follow government regulations and environmental requirements.

EMC AND ELECTRICAL SAFETY EVALUATION: The TRU®/PUR® Surgical Kit does not require EMC and Electrical Safety Evaluation.


MR SAFETY: Sterngold TRU®/PUR® Surgical Kit has not been evaluated for safety and compatibility in the MR environment. The TRU®/PUR® Surgical Kit has not been tested for heating or migration, or image artifact in the MR environment. The safety of TRU®/PUR® Surgical Kit in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

WARRANTY: Sterngold warrants its products to be free from defects in material and/or workmanship. No other warranty is expressed or implied. This warranty applies to the original purchaser only. In the unlikely event of a defect, please follow the returned goods policy outlined in the product catalog.

SYMBOLS: The following table describes the symbols that may be printed on the packaging label. Please refer to the packaging label for the applicable symbols related to the product.

| Symbol | Used For | Symbol | Used For |
|---|------------------------------|--|----------------------------------|
|  | Use by Prescription Only |  | Part No. |
|  | Consult Instructions for Use |  | Batch Code (lot #) |
|  | Do Not re-use |  | Do not use if package is damaged |
|  | Non-Sterile |  | Medical Device |

STERILIZATION: Individual parts should be placed in an appropriate autoclave for moist heat sterilization or in a dry heat pouch for dry heat sterilization. The following sterilization parameters (method, time, and temperature) are required to achieve a 10-6 sterility

| Symbol | Used For | Symbol | Used For |
|---|-------------------------------|---|---|
|  | Sterilized using Irradiation |  | UDI Unique Device Identifier |
|  | Double sterile barrier system |  | Caution, consult accompanying documents |
|  | Do not re-sterilize |  | Use by Date |

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Sterngold Dental, LLC