

**Safety Data Sheet**  
**according to decree (EG) No. 1272 / 2008**

Revision: 19

Revision Date: 31.01.2023, SW

---

**1. Identification of the substance/mixture and of the company/undertaking**

**1.1 Product name (product identifiers):**

**SternTek, Rose, Custom Tray Material**

**1.2 Relevant identified uses of the substance or mixture and uses advised against:**

**Intended purpose:**

Making of individual impression- and functional trays

Occlusal registration

Precized impression for part- or totalprosthetics

**Disadvised uses:**

Impressions

Prosthesis Plastic

Filling Material

**1.3 Details of the Supplier of the safety data sheet:**

**Manufacturer/Company:**

**dent a pharm** Produktionsgesellschaft mbH

**Street:**

Schusterring 35

**Zip-Code / Town**

D-25355 Barmstedt

**Contact for technical information**

+49 (0) 41 23 / 9225 0

**Phone / Fax / E-Mail**

+49 (0) 4123 / 9225 0 / +49 (0) 4123 / 9225 48 / info@dent-a-pharm.de

**1.4 Emergency number**

**Emergency information (24 hours)**

Dr. Wolfgang Willmann

Phone: +49 (0) 170 – 405 30 52

**Safety Data Sheet**  
**according to decree (EG) No. 1272 / 2008**

Revision: 19

Revision Date: 31.01.2023, SW

**2. Hazards Identification****2.1 Classification of the substance or mixture:****According to Regulation (EC) No. 1272/2008:**

- Skin sens. 1 (may cause allergic reactions)

**2.2 Label elements****Labelling according to Regulation 1272/2008 (GHS/CLP)**

Warning

|             |   |
|-------------|---|
| <b>H317</b> | <b>may cause an allergic skin reaction</b>  |
| P272        | contaminated work clothing should not wear outside the workplace  |
| P280        | wear protective gloves/protective clothing/ eye protection/ face protection                                     |
| P302        | in contact with skin:   |
| P352        | flush with plenty of water and soap   |
| P333        | when skin irritation or rash occurs:  |
| P313        | get medical advice / attention  |
| P362        | take off contaminated clothing and wash before wearing again  |
| P364        | wash before wearing again   |
| P501        | ingredients/container have to be disposed in accordance with local/regional/national/ international regulations |

**2.3 Additional hazards:**

- Can become warm during light curing process
- Do not use if known allergy against Acrylates and Methacrylates

**Safety Data Sheet**  
according to decree (EG) No. 1272 / 2008

Revision: 19

Revision Date: 31.01.2023, SW

**3. Composition/Information on Ingredients****Chemical characterization:****Hazardous ingredients:**

| <b>Ingredients</b>                    | <b>CAS-No./ EINECS-No./ EG-No.</b>         | <b>Percentage (%)</b> | <b>Identification</b>   |
|---------------------------------------|--|-----------------------|---|
| Triethylenglykoldimethacrylate        | CAS-No. 109-16-0<br>EINECS-No.203-652-6    | 2-5                   | Skin Sens. 1B H317<br>Skin sensitisation  |
| Phosphine oxide                       | CAS-No. 75980-60-8<br>EG-No. 278-355-8     | 0,05 – 0,1            | Repr 2 H361f<br>suspected human carcinogen; probably harmful to fertility, probably damaging effect on the unborn child<br>Aqu. Chron. 3, H412<br>hazardous to the aquatic environment, chronic |
| Vinylester Resin                      | CAS Nr. 55818-57-0                         | 5-8                   | Skin Sens. 1B H317<br>Skin sensitisation  |
| 2,2-Dimethoxy-2-diphenylacetophenon   | CAS-Nr.24650-42-8<br>EINECS-Nr. 246-386-6  | < 0,1 %               | Aquatic Chronic 1, H410   |
| Campherchinon                         | CAS-Nr. 10373-78-1                         | < 0,1 %               | Skin Irrit. 2, H315<br>STOT SE, H335  |
| Dimethylhydroxyacetophenon            | CAS-Nr. 7473-98-5<br>EINECS-Nr. 231-272-0  | < 0,1 %               | Skin Sens. 1, H317<br>Acute Tox. 4<br>Aquatic Acute   |
| 2-Ethylhexyl-4-(dimethylamino)benzoat | CAS-Nr. 21245-02-3<br>EINECS-Nr. 244-289-3 | < 0,1 %               | Eye Irrit. 2, H319<br>Aquatic Acute, H400<br>Repr. 1B, H360   |
| 2,2-Diethoxyacetophenon               | CAS-Nr. 6175-45-7<br>EINECS-Nr. 228-220-4  | <0,1 %                | STOT SE 3, H335   |

**Safety Data Sheet**  
**according to decree (EG) No. 1272 / 2008**

Revision: 19

Revision Date: 31.01.2023, SW

---

#### **4. First-Aid-Measures**

##### **4.1 Description of first aid measures**

###### **If inhaled:**

Take person to fresh air.  
No further first aid measures are required.

###### **In case of skin contact:**

Wash off skin with soap and water and flush with plenty of Water.  
Contact a physician if the skin is irritated.

###### **In case of eye contact:**

Rinse open eyes an efficient time with plenty of Water.  
Contact ophthalmologist if irritation is persisting.

###### **If swallowed:**

Seek medical advice immediately.

##### **4.2 Most important symptoms and effects, both acute and delayed:**

Sensitization of skin is possible

##### **4.3 Indication of any immediate medical attention and special treatment needed**

Contact physician if sensitization of skin is persisting

#### **5. Fire-Fighting-Measures**

##### **5.1 Extinguishing media**

Suitable extinguishing media:  
Foam, Dry powder, carbon dioxide, waterspray-stream.  
Adjust fire-fighting-measures to ambiance.

For safety reasons unsuitable extinguishing media: none

##### **5.2 Special hazards arising from the substance or mixture**

Particular danger through the product itself, its combustion products or gases are not known

##### **5.3 Advice for fire-fighters**

Use conventional measures to combat chemical fires, that means use self-contained breathing apparatus (isolating-unit).

**Safety Data Sheet**  
**according to decree (EG) No. 1272 / 2008**

Revision: 19

Revision Date: 31.01.2023, SW

---

**6. Accidental release measures**

**6.1 Personal precautions, protective equipment and emergency procedures:**

Arrange for fresh air.

**6.2 Environmental protection measures:**

Do not drain into canalization.

**6.3 Methods for cleaning / absorption:**

Pick up material mechanically. Dispose only complete light cured material. Dispose material under adherence to the local official rules for hazardous wastes.

**6.4 Reference to other sections**

None

**7. Handling and Storage**

**7.1 Precautions for safe handling:**

**Recommendation:**

The conventional precautions for works with chemicals are to be adhered.  
Wear overall while working.

**Information about fire- and explosion protection:**

No special measures necessary.

**7.2 Conditions for safe storage, including any incompatibilities:**

Store at a cool, dry and light protected place.  
Storage time: 3 years if storage conditions will be kept

**7.3 Specific end use(s)**

Wear gloves while working.  
Good ventilation during grinding the cured material is recommended

**Safety Data Sheet**  
**according to decree (EG) No. 1272 / 2008**

Revision: 19

Revision Date: 31.01.2023, SW

## 8. Exposure Controls/Personal Protection

### 8.1 Control parameters

None known

PNCE and DNCL-Values as to be possibly worked out in a preparation safety report do not exist.

### 8.2 Exposure controls

Personal protective equipment::

Pay attention of protection- and hygienic measures. Protection clothes must be stored separately. The usual precaution measures must be attended.

Breathing protection:

not necessary

Eye protection:

Wear protection goggles when working at material

Hand protection:

Sensibilisation: Wearing of applicable protection gloves

is recommended. (by experience of the past/literature) for example Ethiparat®-gloves for 5-10 minutes, application of No-touch-Technik, minimal-contact-techniques is recommended

Body protection:

Wear business usual protection clothes.

Other protection measures:

none

Advices for hygienic  
at work station:

Comply business usual protection- and hygienic  
measures.

Preventive skin protection.

Wash hands and skin after finishing work.

**Safety Data Sheet**  
**according to decree (EG) No. 1272 / 2008**

Revision: 19

Revision Date: 31.01.2023, SW

**9. Physical and chemical Properties****9.1 Information on basic physical and chemical properties**

All Information are for the Acrylate Resin, except point the appearance

|  |  |
|--|--|
| <b>Appearance:</b>                                   | Form: upper jaw and lower jaw pattern, kneadable mass                                |
| <b>Odour/Flavour:</b>                                | acrylates / methacrylates  |
| <b>Odour Threshold:</b>                              | no data available  |
| <b>pH-level:</b>                                     | no data available  |
| <b>Melting point/freezing point:</b>                 | no data available  |
| <b>Initial boiling point and boiling range:</b>      | no data available  |
| <b>Flash point:</b>                                  | not available (approx. > 100°C)<br>appreciated by the flash point of the ingredients |
| <b>Evaporation rate</b>                              | no data available  |
| <b>Flammability:</b>                                 | no data available  |
| <b>Upper/lower flammability or explosive limits:</b> | no data available  |
| <b>Vapour pressure</b>                               | no data available  |
| <b>Vapour density:</b>                               | no data available  |
| <b>Relative density:</b>                             | no data available  |
| <b>Water solubility:</b>                             | unsoluble in water   |
| <b>Partition coefficient:</b>                        |  |
| <b>n-Octanol/water:</b>                              | no data available  |
| <b>Auto ignition temperature:</b>                    | no data available  |
| <b>Decomposition temperature:</b>                    | no data available  |
| <b>Viscosity:</b>                                    | no data available  |
| <b>Explosive properties</b>                          | no data available  |
| <b>Oxidizing properties</b>                          | no data available  |
| <b>9.2 Other safety information:</b>                 | no data available  |

**10. Stability and reactivity****10.1 Reactivity**

Product polymerizes slowly by contact with Sunlight or other Light

**10.2 Chemical stability**

If protected against heat, solar radiation and light the material is chemical stable.

**10.3 Possibility of hazardous reactions**

No hazardous heat emission during slow polymerizing

**10.4 Conditions to avoid**

Protect against Sunlight, heat and other light to avoid slowly polymerization.

**10.5 Incompatible materials**

Protect against radical formers to avoid slowly polymerization

**10.6 Hazardous decomposition products:**

None, if used according to directions.

After considerably exceedance of storage time and / or exceedance of storage temperature it is possible that the product starts polymerizing under low heat impact.

**11. Toxicological information**

**Safety Data Sheet**  
**according to decree (EG) No. 1272 / 2008**

Revision: 19

Revision Date: 31.01.2023, SW

---

**11.1 Information of toxicological effects**

**Acute Toxicity**

According to the single ingredients the mixture is not acute toxic

**Skin corrosion/irritation**

No irritation to skin

**Serious eye damage/eye irritation**

No irritation to eyes

**Respiratory or skin irritation**

No irritation to Respiratory or skin

**Skin sensitisation**

Skin sensitisation possible according to the ingredients

**Toxicity – repeated exposure**

No toxicity known after repeated exposure

**Germ cell mutagenicity**

**Carcinogenicity**

No carcinogenicity according to the ingredients

**Mutagenicity**

No mutagenicity according to the ingredients

**Reproductive toxicity**

No reproductive toxicity according to the concentration of the single components.

**Kinetic of the toxic reaction**

According the new Safety Data Sheet according to REACH, specifications about kinetic of the toxic reaction, metabolism of the substance and allocation of the substance in the human body must be mentioned at this place. Data for this is not available .

**Further Information**

By appropriate exposure and conventionally use and to our experiences and informations the product causes no hazardous effects.

The product has to handled with the usual care when handling chemicals.

Methacrylates may cause allergical reactions.



**Safety Data Sheet**  
**according to decree (EG) No. 1272 / 2008**

Revision: 19

Revision Date: 31.01.2023, SW

**12. Ecological information****12.1 Toxicity**

Fish toxicity: no data available for this mixture

Water organism: no data available

Bacteria toxicity: no data available

**12.2 Persistence and degradability**

The ingredients of the mixture are not degradable

**12.3 Bioaccumulation potential:**

No data available

**12.4 Mobility in soil**

No data available

**12.5 Effect of the determination of the PTB- and vPvB-Report:**

If the material causes hormonal effects or prevents them is unknown to us.

**12.6 Other adverse effects**

No data available

**13. Disposal considerations****13.1 Waste treatment methods**

To dispose of accordingly with notice of the local regulations of authority.  
Do not empty into drains.

**Disposal key according to AVV regulation**

16 03 05 (organic disposal, that contains hazardous substances)

**Package**

Packages that are not contaminated can be recycled.

**14. Transport information/rules****14.1 UN number**

none

**14.2 UN proper shipping name**

none

**14.3 Transport hazard class(es)**

none

**14.4 Packing group**

none

**14.5 Environmental hazards**

none

**14.6 Special precautions for Transport**

none

**14.7 Transportation in bulk according to Annex II of MARPOL 73/78 and IBC Code**

no

**Safety Data Sheet**  
**according to decree (EG) No. 1272 / 2008**

Revision: 19

Revision Date: 31.01.2023, SW

---

**15. Regulatory Information****15.1 Safety health and environmental/legislation specific for the substances or mixture****GHS-Labeling according to guideline 1272/2008 (GHS/CL)**

Warning

|             |   |
|-------------|---|
| <b>H317</b> | <b>may cause an allergic skin reaction</b>  |
| P272        | contaminated work clothing should not wear outside the workplace  |
| P280        | wear protective gloves/protective clothing/eye protection/<br>face protection                                 |
| P302        | in contact with skin:   |
| P352        | flush with plenty of water and soap   |
| P333        | when skin irritation or rash occurs:  |
| P313        | get medical advice / attention  |
| P362        | take off contaminated clothing and wash before wearing again  |
| P364        | wash before wearing again   |
| P501        | ingredients/container have to be disposed in accordance with local/regional/national/international regulation |

**15.2 Safety Assessments**

Safety Assessments for substances in this preparation were not executed

**Safety Data Sheet**  
**according to decree (EG) No. 1272 / 2008**

Revision: 19

Revision Date: 31.01.2023, SW

**16. Other information****Further applicable EG Guidelines**

Guideline for preparation (1999/45/EC), at last changed by the guideline 2006/8/EC.  
Substance Guideline (67/548/EWG), at last changed by the guideline 2006/121 EC.  
REACH-Guideline (EC) No. 1907/2006, at last changed by the guideline 1272/2008 EC.

**Manufacturer's suggested application restriction:**

Only for commercial use – no public product

**H-phrases on which is referred in chapter 2 and 3**

H317            May cause an allergic skin reaction  
H361f          May impair fertility  
H412            Harmful to aquatic life with long lasting effects.

**Changes compared to the last version**

Change of Labelling according GHS-Labelling

**Department issuing SDS:**

Dr. Wolfgang Willmann (Phone: +49 (0) 41 23 – 9225 0)

**Additional Information:**Origin: <sup>1</sup> <http://www.baua.de><sup>2</sup> <http://www.arbeitssicherheit.de>

Some of the information presented and conclusions drawn herein are from sources other than direct test-data on the product itself.

The contents and format of this SDS are in accordance with EEC Commission (EC) No. 1272 / 2008.

**Disclaimer of Liability**

The information in this SDS was obtained from sources which we believe are reliable. However, the information is provided without any warranty, express or implied, regarding its correctness.

The conditions or methods of handling, storage, use or disposal of the product are beyond our control and may be beyond our knowledge.

For this and other reasons, we do not assume responsibility and expressly disclaim liability for loss, damage or expense arising out of or in any way connected with the handling, storage, use or disposal of the product.

This SDS was prepared and is to be used only for this product. If the product is used as a component in another product, this SDS information may not be applicable.