



23 Frank Mossberg Drive · Attleboro, MA 02703

Tel: (508) 226-5660
Cust. Serv: (800) 243-9942
Fax: (508) 226-5473
Toll Free Fax: (800) 531-2685

Alloys
Attachments
Implants <http://www.sterngold.com>
Restorative
Systems

MATERIAL SAFETY DATA SHEET

Sterngold Part Number: 811165, 811017, 811360, 811039

Sterngold Description: ERA nylon grey males

Supplier Number: AB72620234

Supplier Description: NY Zytel ST801 025.00:1 Grey

Supplier MSDS Identification: n/a

1/24/12

Date Updated:

M A T E R I A L S A F E T Y D A T A S H E E T

DATE: 11/09/11
PAGE: 1

SHIPPED FROM: CLARIANT CORPORATION, MASTERBATCHES DIVISION
HOLDEN-CLARIANT MASTERBATCHES
85 HOLDEN INDUSTRIAL PARK
HOLDEN MA 01520
PHONE: 508-829-6321

EMERGENCY TELEPHONE: 269-753-8434

CLARIANT WEBSITE: HTTP:\\WWW.CLARIANTMASTERBATCH.COM

PRODUCT DESCRIPTION: NY ZYTEL ST801 025.00:1 GRAY
PRODUCT CODE: AB72620234

***** MATERIAL SAFETY DATA SHEET *****

*
THE ATTACHED MATERIAL SAFETY DATA SHEET (MSDS) IS AN IMPORTANT DOCUMENT.
THE OSHA HAZARDOUS COMMUNICATION STANDARD REQUIRES US TO FORWARD A MSDS
PRIOR TO OR ACCOMPANYING AN INITIAL SHIPMENT OF A PRODUCT AND THE
"STANDARD" REQUIRES YOUR COMPANY TO KEEP THE MSDS ON FILE.
*

*
IF YOU ARE NOT RESPONSIBLE FOR HEALTH AND SAFETY IN YOUR COMPANY PLEASE
INSURE THAT YOUR HEALTH AND SAFETY MANAGER RECEIVES THIS DOCUMENT.
*

*
THANK YOU.

HAZARDOUS MATERIAL IDENTIFICATION SYSTEM (HMIS)

HEALTH ----->1
FLAMMABILITY ----->1
REACTIVITY ----->0
PERSONAL PROTECTION -->B

M A T E R I A L S A F E T Y D A T A S H E E T

DATE: 11/09/11
PAGE: 2

SHIPPED FROM: CLARIANT CORPORATION, MASTERBATCHES DIVISION
HOLDEN-CLARIANT MASTERBATCHES
85 HOLDEN INDUSTRIAL PARK
HOLDEN MA 01520
PHONE: 508-829-6321

EMERGENCY TELEPHONE: 269-753-8434

CLARIANT WEBSITE: HTTP:\\WWW.CLARIANTMASTERBATCH.COM

PRODUCT DESCRIPTION: NY ZYTEL ST801 025.00:1 GRAY
PRODUCT CODE: AB72620234

SECTION I - PRODUCT INFORMATION:

PRODUCT DESCRIPTION: NY ZYTEL ST801 025.00:1 GRAY
PRODUCT CODE: AB72620234 VERSION NUMBER: 10 DATE: 11/09/11
CHEMICAL NAME/FAMILY:
THERMOPLASTICS RESIN(S) + COLORANT

SECTION II - HAZARDOUS INGREDIENTS/IDENTITY INFORMATION:

THIS PRODUCT CONTAINS, AS A COMPONENT, THE SUBSTANCE(S) LISTED BELOW
WHICH IS IDENTIFIED AS A HAZARDOUS CHEMICAL UNDER THE CRITERIA OF THE
OSHA HAZARD COMMUNICATION STANDARD 29CFR1910.1200

CAS NUMBER	DESCRIPTION	OSHA		OTHER LIMITS	DOES NOT EXCEED (%)
		PEL MG/M3	TLV MG/M3		
68186-90-3	CHROME ANTIMONY TITANIUM	-	-	-	1.23
CRCPD	CHROMIUM (III) COMPOUND	0.5	0.5	-	----
SBCMPD	ANTIMONY & COMPOUNDS (AS SB)	0.5	0.5	-	----
1345-16-0	COBALT ALUMINATE	-	-	-	.23
68186-91-4	COPPER CHROMITE SPINEL	N.E.	N.E.	-	6.88
13463-67-7	TITANIUM DIOXIDE	15.0	10.0	-	28.64
1344-28-1	ALUMINUM OXIDE (DUST)	10.0	10.0	5.0	.29

SECTION III - PHYSICAL/CHEMICAL CHARACTERISTICS:

BOILING POINT: NA
SPECIFIC GRAVITY: 1.13
VAPOR PRESSURE: NA
MELTING POINT: Non Determined
VAPOR DENSITY: NA
EVAPORATION RATE: NA
SOLUBILITY IN WATER: INSOLUBLE
APPEARANCE AND ODOR: COLORED PLASTIC PELLETT - NO SIGNIFICANT ODOR

SECTION IV - FIRE AND EXPLOSION HAZARD DATA:

FLASH POINT (METHOD USED):
THERE IS A COMPONENT IN THIS PRODUCT WITH A FLASH POINT OF 350 C.
THE METHOD USED TO DETERMINE THIS IS: NA

M A T E R I A L S A F E T Y D A T A S H E E T

DATE: 11/09/11
PAGE: 3

SHIPPED FROM: CLARIANT CORPORATION, MASTERBATCHES DIVISION
HOLDEN-CLARIANT MASTERBATCHES
85 HOLDEN INDUSTRIAL PARK
HOLDEN MA 01520
PHONE: 508-829-6321

EMERGENCY TELEPHONE: 269-753-8434

CLARIANT WEBSITE: HTTP:\\WWW.CLARIANTMASTERBATCH.COM

PRODUCT DESCRIPTION: NY ZYTEL ST801 025.00:1 GRAY
PRODUCT CODE: AB72620234

SECTION IV - FIRE AND EXPLOSION HAZARD DATA: (CONT)

ALL OTHER COMPONENTS USED IN THIS PRODUCT EITHER HAVE A HIGHER FLASH
POINT OR ARE NOT COMBUSTIBLE AND DO NOT HAVE A FLASH POINT.

EXTINGUISHING MEDIA:

WATER SPRAY, FOAM, CARBON DIOXIDE, OR DRY CHEMICAL

SPECIAL FIRE FIGHTING PROCEDURES:

FIREFIGHTERS SHOULD BE EQUIPPED WITH SELF-CONTAINED BREATHING APPARATUS TO
PROTECT AGAINST POTENTIALLY TOXIC AND IRRITATING FUMES. AVOID DUSTING.

UNUSUAL FIRE AND EXPLOSION HAZARDS:

COMPLETE COMBUSTION IS EXPECTED TO PRODUCE CARBON, OXIDES OF CARBON AND
NITROGEN, PLUS TRACE AMOUNTS OF ELEMENTAL CHEMICALS.

SECTION V - REACTIVITY DATA:

STABILITY: STABLE -> X UNSTABLE ->

STABILITY (CONDITIONS TO AVOID):

AVOID EXPOSURE TO HEAT, SPARK, OR OPEN FLAME.

INCOMPATIBILITY (MATERIALS TO AVOID):

AVOID CONTACT WITH STRONG ALKALI, AMINES.

HAZARDOUS DECOMPOSITION OR BYPRODUCTS:

THERMAL DECOMPOSITION PRODUCTS MAY INCLUDE HAZARDOUS FUMES OF ANTIMONY
OXIDE, MANGANESE, AND TITANIUM DIOXIDE.

HAZARDOUS POLYMERIZATION: WILL OCCUR -> WILL NOT OCCUR -> X

HAZARDOUS POLYMERIZATION (CONDITIONS TO AVOID):

AVOID CONTACT WITH STRONG CHEMICAL OXIDIZING AGENTS.

SECTION VI - HEALTH HAZARD DATA:

HEALTH HAZARDS (ACUTE AND CHRONIC):

HEALTH HAZARD: INHALATION

M A T E R I A L S A F E T Y D A T A S H E E T

DATE: 11/09/11

PAGE: 4

SHIPPED FROM: CLARIANT CORPORATION, MASTERBATCHES DIVISION
HOLDEN-CLARIANT MASTERBATCHES
85 HOLDEN INDUSTRIAL PARK
HOLDEN MA 01520
PHONE: 508-829-6321

EMERGENCY TELEPHONE: 269-753-8434

CLARIANT WEBSITE: HTTP:\\WWW.CLARIANTMASTERBATCH.COM

PRODUCT DESCRIPTION: NY ZYTEL ST801 025.00:1 GRAY
PRODUCT CODE: AB72620234

SECTION VI - HEALTH HAZARD DATA: (CONT)

SHORT TERM OVEREXPOSURE BY IHALATION TO TITANIUM DIOXIDE MAY CAUSE IRRITATION OF NOSE, THROAT, AND LUNGS WITH COUGH, DIFFICULTY BREATHING OR SHORTNESS OF BREATH. STUDIES SHOWED THAT EMPLOYEES WHO HAD BEEN EXPOSED TO TITANIUM DIOXIDE WERE AT NO GREATER RISK OF DEVELOPING LUNG CANCER THAN WERE EMPLOYEES WHO HAD NOT BEEN EXPOSED TO TITANIUM DIOXIDE. INHALATION OF AMORPHOUS SILICA MAY CAUSE DRYING OF MUCOUS MEMBRANES AND IRRITATION OF NOSE, THROAT, AND LUNGS WITH NOSEBLEEDS, COUGH AND DIFFICULTY BREATHING OR SHORTNESS OF BREATH.

HEALTH HAZARD: SKIN
MAY CAUSE IRRITATION TO SKIN

HEALTH HAZARD: INGESTION
EXPECTED TO BE A LOW INGESTION HAZARD

HEALTH HAZARD: EYE CONTACT
CONTACT WITH EYE MAY CAUSE EYE IRRITATION

IS ANY CHEMICAL IN THIS PRODUCT AT A CONCENTRATION OF 0.1 % OR MORE CLASSIFIED AS A CARCINOGEN BY IARC, NTP, OR OSHA? YES

SIGNS AND SYMPTOMS OF EXPOSURE:

* THE FOLLOWING INFORMATION IS PROVIDED ON THE BELOW HAZARDOUS INGREDIENT(S) PER OSHA HAZARD COMMUNICATION STANDARD 29CFR 1910.1200.

CHROME ANTIMONY TITANIUM BUFF RUTILE: REPEATED OVEREXPOSURE MAY CAUSE EYE, SKIN, AND RESPIRATORY TRACT IRRITATION. SOME OF THE COMPOUNDS HAVE DEMONSTRATED VARIOUS TOXIC PROPERTIES. ROUTES OF ENTRY: INHALATION AND INGESTION.

COBALT ALUMINATE: IRRITATION OF THE EYES. SKIN AND RESPIRATORY.
IARC HAS CLASSIFIED COBALT AND COBALT COMPOUNDS AS POSSIBLY CARCINOGENIC TO HUMANS (GROUP 2B)

COPPER CHROMITE: OVEREXPOSURE MAY CAUSE EYE, SKIN AND RESPIRATORY TRACT IRRITATION.

M A T E R I A L S A F E T Y D A T A S H E E T

DATE: 11/09/11
PAGE: 5

SHIPPED FROM: CLARIANT CORPORATION, MASTERBATCHES DIVISION
HOLDEN-CLARIANT MASTERBATCHES
85 HOLDEN INDUSTRIAL PARK
HOLDEN MA 01520
PHONE: 508-829-6321

EMERGENCY TELEPHONE: 269-753-8434

CLARIANT WEBSITE: HTTP:\\WWW.CLARIANTMASTERBATCH.COM

PRODUCT DESCRIPTION: NY ZYTEL ST801 025.00:1 GRAY
PRODUCT CODE: AB72620234

SECTION VI - HEALTH HAZARD DATA: (CONT)

TITANIUM DIOXIDE, ALUMINUM OXIDE: TITANIUM DIOXIDE HAS BEEN LISTED BY IARC AS A GROUP 2(B) POSSIBLE CARCINOGEN TO HUMANS. MAY CAUSE IRRITATION TO EYES, SKIN AND RESPIRATORY TRACT AND MAY AGGRAVATE PRE-EXISTING RESPIRATORY CONDITIONS. MAY ALSO CAUSE TEMPORARY DRYING EFFECT.

MEDICAL CONDITIONS GENERALLY AGGRAVATED BY EXPOSURE:

* THE FOLLOWING INFORMATION PROVIDED ON THE BELOW HAZARDOUS INGREDIENT(S) APPLIES TO THE INGREDIENT(S) IN THE PURE DRY FORM. THE FORM OF THE INGREDIENT(S) PROVIDED TO YOU IS ENCAPSULATED IN PLASTIC, AND THEREFORE, THE LIKELIHOOD OF EXPOSURE IS MUCH LESS, EVEN TO THE POINT OF NEGLIGIBLE. WE THEREFORE CONSIDER THE LEVEL OF HAZARD TO BE NEGLIGIBLE.

CHROME ANTIMONY TITANIUM BUFF RUTILE: RESPIRATORY AND SKIN DISORDERS.

COBALT ALUMINATE: RESPIRATORY AND SKIN DISORDERS AGGRAVATED BY DUST.

COPPER CHROMITE: RESPIRATORY AND SKIN DISORDERS AGGRAVATED BY DUST.

TITANIUM DIOXIDE AND ALUMINUM OXIDE: MAY AGGRAVATE PRE-EXISTING RESPIRATORY CONDITIONS.

EMERGENCY FIRST AID: INHALATION

REMOVE TO UNCONTAMINATED AREA. IF NOT BREATHING, ADMINISTER ARTIFICIAL RESPIRATION. IF BREATHING IS DIFFICULT, GIVE OXYGEN. SEEK MEDICAL ATTENTION

EMERGENCY FIRST AID: SKIN

REMOVE CONTAMINATED CLOTHING AND SHOES IMMEDIATELY. WASH WITH SOAP OR MILD DETERGENT AND LARGE AMOUNTS OF WATER UNTIL NO EVIDENCE OF CHEMICAL REMAINS (AT LEAST 15 - 20 MINUTES). GET MEDICAL ATTENTION IMMEDIATELY.

EMERGENCY FIRST AID: INGESTION

FOR EXCESSIVE INGESTION DRINK WATER OR MILK. CALL A PHYSICIAN

EMERGENCY FIRST AID: EYE CONTACT

WASH EYES IMMEDIATELY WITH LARGE AMOUNTS OF WATER OR NORMAL SALINE, OCCASIONALLY LIFTING UPPER AND LOWER LIDS, UNTIL NO EVIDENCE OF CHEMICAL REMAINS (AT LEAST 15 - 20 MINUTES). GET MEDICAL ATTENTION IMMEDIATELY.

M A T E R I A L S A F E T Y D A T A S H E E T

DATE: 11/09/11
PAGE: 6

SHIPPED FROM: CLARIANT CORPORATION, MASTERBATCHES DIVISION
HOLDEN-CLARIANT MASTERBATCHES
85 HOLDEN INDUSTRIAL PARK
HOLDEN MA 01520
PHONE: 508-829-6321

EMERGENCY TELEPHONE: 269-753-8434

CLARIANT WEBSITE: HTTP:\\WWW.CLARIANTMASTERBATCH.COM

PRODUCT DESCRIPTION: NY ZYTEL ST801 025.00:1 GRAY
PRODUCT CODE: AB72620234

SECTION VII - PRECAUTIONS FOR SAFE HANDLING AND USE:

HANDLING: SPILLAGE

SWEEP UP AND PLACE IN A SUITABLE CONTAINER FOR RECLAMATION OR DISPOSAL.

HANDLING: WASTE DISPOSAL

LANDFILL OR INCINERATE IN ACCORDANCE WITH LOCAL, STATE AND FEDERAL LAWS.

HANDLING: STORAGE

AVOID GENERATING DUST. AVOID CONTACT WITH EYES, SKIN, OR CLOTHING. WASH THOROUGHLY AFTER HANDLING. STORE IN A DRY WELL VENTILATED AREA. KEEP CONTAINER TIGHTLY CLOSED.

HANDLING: OTHER PRECAUTIONS

STORE MATERIAL IN A CLOSED CONTAINER. NORMAL WAREHOUSING.

SECTION VIII - CONTROL MEASURES:

CONTROL: RESPIRATORY PROTECTION

UNLESS AIR MONITORING DEMONSTRATES DUSTS/MIST/VAPORS ARE BELOW THE RECOMMENDED EXPOSURE LEVELS, WEAR A PROPERLY FITTED RESPIRATOR (NIOSH/MSHA) OR DUST MASK DURING EXPOSURE.

CONTROL: VENTILATION

USE LOCAL VENTILATION

CONTROL: GLOVES

IMPERVIOUS GLOVES

CONTROL: EYE PROTECTION

GOGGLES

CONTROL: OTHER CLOTHING

TO PREVENT SKIN CONTACT WEAR SKIN PROTECTION, SUCH AS IMPERVIOUS GLOVES, APRON, WORKPANTS, LONG SLEEVE WORKSHIRT, OR DISPOSABLE COVERALLS.

CONTROL: WORK/HYGIENE

OBSERVE PROPER INDUSTRIAL HYGIENIC PRACTICES ACCORDING TO REGULATIONS.

SECTION IX - OTHER REGULATORY INFORMATION:

M A T E R I A L S A F E T Y D A T A S H E E T

DATE: 11/09/11
PAGE: 7

SHIPPED FROM: CLARIANT CORPORATION, MASTERBATCHES DIVISION
HOLDEN-CLARIANT MASTERBATCHES
85 HOLDEN INDUSTRIAL PARK
HOLDEN MA 01520
PHONE: 508-829-6321

EMERGENCY TELEPHONE: 269-753-8434

CLARIANT WEBSITE: HTTP:\\WWW.CLARIANTMASTERBATCH.COM

PRODUCT DESCRIPTION: NY ZYTEL ST801 025.00:1 GRAY
PRODUCT CODE: AB72620234

SECTION IX - OTHER REGULATORY INFORMATION:

SARA - THIS PRODUCT CONTAINS THE FOLLOWING TOXIC CHEMICAL(S) SUBJECT TO THE REPORTING OF SECTION 313 OF TITLE III OF THE SUPERFUND AMENDMENTS AND REAUTHORIZATION ACT OF 1986 AND 40 CFR PART 372: THIS INFORMATION SHOULD BE INCLUDED IN ALL MSDS'S THAT ARE COPIED AND DISTRIBUTED FOR THIS MATERIAL.

CAS #	DOES NOT EXCEED (%)	DESCRIPTION
SBCMPD	1.23	ANTIMONY COMPOUND
CRCMPD	8.11	CHROMIUM (III) COMPOUND
7440-47-3	3.02	CHROMIUM (III) AS (CR)
7440-36-0	.15	ANTIMONY (SB) METAL
COCMPD	.23	COBALT COMPOUND
7440-48-4	.06	COBALT METAL (AS CO)
CUCMPD	6.88	COPPER COMPOUND
7440-50-8	2.06	COPPER (DUST)

REACH REGULATION - THE SVHC REPORTED IN THE PROPOSAL LIST FOR INCLUSION IN ANNEX XIV (ART 59.4) AND ACCORDING TO OUR KNOWLEDGE AND THE REACH REGULATION ITSELF. DETAILS ABOUT SVHC CONTAINED IN OUR MB ARE HEREUNDER REPORTED AS OF FEBRUARY 2010.

CAS #	DOES NOT EXCEED (%)	DESCRIPTION
NONE		

TSCA - ALL INGREDIENTS USED TO PRODUCE THE ABOVE FORMULATION ARE IN COMPLIANCE WITH THE TOXIC SUBSTANCES AND CONTROL ACT INVENTORY AS REQUIRED BY FEDERAL LAW.

SECTION X - STATE REGULATORY INFORMATION:

CA PROPOSITION 65: WARNING - ANY COMPONENT OF THIS PRODUCT WHICH CONTAINS MATERIALS KNOWN TO THE STATE OF CALIFORNIA TO CAUSE CANCER ARE LISTED BELOW:

PLEASE NOTE IT IS THE RESPONSIBILITY OF THE END USER TO DETERMINE IF ITS FINISHED PRODUCTS REQUIRE LABELING UNDER CALIFORNIA PROPOSITION 65. IMPURITIES OF HEXAVALENT CHROMIUM AND LEAD.

CA PROPOSITION 65: WARNING - ANY COMPONENT OF THIS PRODUCT WHICH CONTAINS MATERIALS KNOWN TO THE STATE OF CALIFORNIA TO

M A T E R I A L S A F E T Y D A T A S H E E T

DATE: 11/09/11
PAGE: 8

SHIPPED FROM: CLARIANT CORPORATION, MASTERBATCHES DIVISION
HOLDEN-CLARIANT MASTERBATCHES
85 HOLDEN INDUSTRIAL PARK
HOLDEN MA 01520
PHONE: 508-829-6321

EMERGENCY TELEPHONE: 269-753-8434

CLARIANT WEBSITE: HTTP:\\WWW.CLARIANTMASTERBATCH.COM

PRODUCT DESCRIPTION: NY ZYTEL ST801 025.00:1 GRAY
PRODUCT CODE: AB72620234

SECTION X - STATE REGULATORY INFORMATION: (CONT)

CAUSE BIRTH DEFECTS OR OTHER REPRODUCTIVE HARM
ARE LISTED BELOW:

PLEASE NOTE IT IS THE RESPONSIBILITY OF THE END USER TO DETERMINE IF ITS
FINISHED PRODUCTS REQUIRE LABELING UNDER CALIFORNIA PROPOSITION 65.

NONE

CANADIAN WHMIS: ANY COMPONENT OF THIS PRODUCT WHICH IS A HAZARDOUS CHEMICAL
AS GIVEN IN WHMIS IS SHOWN BELOW:

ANTIMONY - D1B - TOXIC MATERIAL CAUSING IMMEDIATE AND SERIOUS TOXIC EFFECTS
CLASS D DIV. 2B
TITANIUM DIOXIDE: D2B

FLORIDA: NON-HAZARDOUS INGREDIENTS ARE BEING WITHHELD AS TRADE SECRET
INFORMATION. ANY COMPONENT OF THIS PRODUCT WHICH IS A
HAZARDOUS CHEMICAL AS LISTED ON THE FLORIDA HAZARDOUS LIST
IS GIVEN IN SECTION II OF THIS MSDS. ANY LISTED CHEMICALS NOT
SHOWN IN SECTION II ARE AS FOLLOWS:

NONE

MASSACHUSETTS:NON-HAZARDOUS INGREDIENTS ARE BEING WITHHELD AS TRADE SECRET
INFORMATION. ANY COMPONENT OF THIS PRODUCT WHICH IS A
HAZARDOUS CHEMICAL AS LISTED ON THE MASSACHUSETTS HAZARDOUS
LIST IS GIVEN IN SECTION II OF THIS MSDS. ANY LISTED
CHEMICALS NOT SHOWN IN SECTION II ARE AS FOLLOWS:

NONE

NEW JERSEY: NON-HAZARDOUS INGREDIENTS ARE BEING WITHHELD AS TRADE SECRET
INFORMATION. ANY COMPONENT OF THIS PRODUCT WHICH IS A
HAZARDOUS CHEMICAL AS LISTED ON THE NEW JERSEY HAZARDOUS LIST
IS GIVEN IN SECTION II OF THIS MSDS. ANY LISTED CHEMICALS NOT
SHOWN IN SECTION II ARE AS FOLLOWS:

NONE

PENNSYLVANIA: NON-HAZARDOUS INGREDIENTS ARE BEING WITHHELD AS TRADE SECRET
INFORMATION. ANY COMPONENT OF THIS PRODUCT WHICH IS A
HAZARDOUS CHEMICAL AS LISTED ON THE PENNSYLVANIA HAZARDOUS
LIST IS GIVEN IN SECTION II OF THIS MSDS. ANY LISTED
CHEMICALS NOT SHOWN IN SECTION II ARE AS FOLLOWS:

NONE

M A T E R I A L S A F E T Y D A T A S H E E T

DATE: 11/09/11
PAGE: 9

SHIPPED FROM: CLARIANT CORPORATION, MASTERBATCHES DIVISION
HOLDEN-CLARIANT MASTERBATCHES
85 HOLDEN INDUSTRIAL PARK
HOLDEN MA 01520
PHONE: 508-829-6321

EMERGENCY TELEPHONE: 269-753-8434

CLARIANT WEBSITE: HTTP:\\WWW.CLARIANTMASTERBATCH.COM

PRODUCT DESCRIPTION: NY ZYTEL ST801 025.00:1 GRAY
PRODUCT CODE: AB72620234

SECTION X - STATE REGULATORY INFORMATION: (CONT)

RHODE ISLAND: ANY COMPONENT OF THIS PRODUCT WHICH IS A HAZARDOUS CHEMICAL
CHEMICAL AS LISTED IN THE RHODE ISLAND SUBSTANCE LIST IS
SHOWN IN SECTION II OF THIS MSDS. ANY LISTED CHEMICALS NOT
SHOWN IN SECTION II ARE AS FOLLOWS:
NONE

THIS INFORMATION IS ACCURATE AND RELIABLE TO THE BEST OF OUR KNOWLEDGE.
IT IS FURNISHED WITHOUT WARRANTY, EXPRESSED OR IMPLIED. CLARIANT
CORPORATION ASSUMES NO LEGAL RESPONSIBILITY FOR THE USE OF, OR RELIANCE
UPON, THESE DATA FOR HAZARDS WHICH MIGHT BE ASSOCIATED WITH THE USE OF THIS
MATERIAL OR FOR THE RESULTS OBTAINED. IT IS THE RESPONSIBILITY OF THE USER
TO COMPLY WITH ALL APPLICABLE LAWS AND REGULATIONS.